

## DATA SHEET FRONTLINE Plus - Bleed Control Trauma Kit

SKU Number **SS-6010125045**  
 EAN13 Barcode **5060660858170**

Date of Publication **09-01-2025**  
 Commodity Code **30065000.00**



|                            |  |
|----------------------------|--|
| <p>Product Description</p> | <p><b>FRONTLINE Plus – Bleed Control Trauma Kit</b></p> <p>When severe haemorrhage control is needed, having a kit equipped with the right medical products is essential. Simple, user-friendly Bleed Control Kits are critical to the effectiveness of first responders.</p> <p>The <b>STEIN FRONTLINE Plus</b> Bleed Control Trauma Kit includes high-quality medical supplies necessary to manage life-threatening haemorrhages.</p> <p>This kit allows trained First Aiders, Advanced Professionals, and First Responders to provide lifesaving care in critical situations whilst waiting for Emergency Services to arrive.</p> <p>Contained in a durable, lightweight bag made from advanced materials, this kit is designed for quick access in emergencies.</p> <p>The <b>STEIN FRONTLINE Plus</b> Bleed Control Trauma Kit functions as a grab bag that can be easily attached to the harness of a climber performing an Aerial Rescue. For added convenience, it can also be clipped near the casualty, allowing rapid access to its contents for immediate treatment.</p> |
|----------------------------|--|

|                 |   |
|-----------------|---|
| <b>Contents</b> | <ul style="list-style-type: none"> <li>1 x SAM XT Tourniquet</li> <li>1 x SWAT-T Rescue Orange Tourniquet</li> <li>1 x 6" FCP Military Emergency Bandage</li> <li>1 x ChitoSAM Haemostatic Gauze 4' Z-Fold</li> <li>1 x Medical Gauze No14 - Wound Dressing</li> <li>1 x Large First Aid Dressing 18cm x 18cm</li> <li>1 x Burns Dressing 10cm x 10cm</li> <li>1 x Emergency Foil Thermal Blanket</li> <li>1 x Sterile Saline Spray (100ml)</li> <li>4 x 20ml Ampoules Sterile Wound Wash</li> <li>1 x Resuscitation Face Shield</li> <li>1 x Pack of 10 Plasters</li> <li>1 x Pack Sterile Skin Closure Strips (5)</li> <li>1 x High Performance Medical Shears</li> <li>1 x Medical Marker Pen</li> <li>2 x Disposable Medical Gloves</li> <li>2 x De-Contamination Wipes 30cm x 20cm</li> <li>2 x Sterile Wound Cleansing Wipes</li> <li>3 x Security Seals</li> </ul> |
|-----------------|---|

|                         |   |
|-------------------------|---|
| <b>Safety Standards</b> | <ul style="list-style-type: none"> <li>• All medical products used in the STEIN Kits are CE Marked and conform with the relevant European Union harmonisation legislation.</li> </ul> |
|-------------------------|---|

|                       |  |
|-----------------------|--|
| <b>Product Weight</b> | <ul style="list-style-type: none"> <li>• 1125g</li> </ul>                          |
| <b>Size</b>           | <ul style="list-style-type: none"> <li>• 270cm x 170cm x 90cm (approx.)</li> </ul> |

|                 |  |
|-----------------|--|
| <b>Storage</b>  | <ul style="list-style-type: none"> <li>• Keep stored and transported at ambient temperature (0°C to + 30°C)</li> <li>• Do not store in direct/high heat or sunlight as this may affect the internal components of the kit</li> </ul>   |
| <b>Lifetime</b> | <ul style="list-style-type: none"> <li>• Not all products have an expiry date, however most sterile items do. If a non-sterile product has no expiry date and still sealed in its original packaging, it's a question of judgement as to whether it is still suitable for use.</li> <li>• Items that have exceeded their expiry date should be disposed and replaced as they may not be effective and can start to deteriorate, making them unreliable.</li> </ul> |
| <b>Cleaning</b> | <ul style="list-style-type: none"> <li>• Clean the storage bag using warm soapy water and soft cloth</li> </ul>  |
| <b>Disposal</b> | <ul style="list-style-type: none"> <li>• Recycling is preferred where possible. Look out for the recycling symbol for material categories</li> </ul>   |

|                               |   |
|-------------------------------|---|
| <b>Replacement Components</b> | <ul style="list-style-type: none"><li>• Following the use of any product the kit should be replenished with items used so they are fully functional for future incidents</li><li>• Refer to the STEIN Website for replacement items <a href="http://www.steinworldwide.com">www.steinworldwide.com</a></li></ul>                |
| <b>Additional Information</b> | <ul style="list-style-type: none"><li>• In extreme cases of haemorrhage further products may be required over and above the items in these kits.</li><li>• This kit does not replace the mandatory Workplace First Aid Kits complying with BS 8599-1:2019</li></ul>   |
| <b>Training</b>               | <ul style="list-style-type: none"><li>• It is not feasible to address every potential scenario involving the use of this product. Users are strongly advised to obtain professional training from a fully qualified and competent instructor before use.</li><li>• Training Providers can be recommended upon request</li></ul> |
| <b>Trademarks</b>             | <ul style="list-style-type: none"><li>• STEIN® is a registered trademark</li></ul>  |